

# baebies

## SEEKER™

### INSTRUMENT SPECIFICATIONS

Throughput per Day*	960 Samples
Punch-to-Result	3 hours 40 minutes (hands on time = 20 minutes)
Size - Lid Open (W x H x D)	7.75in x 21.75in x 20.75in (20cm x 55cm x 53cm)
Size - Lid Closed (W x H x D)	7.75in x 12in x 20.75in (20cm x 30cm x 53cm)
Weight	23lb (10.5kg)
Detector	Fluorometric
Input Voltage	100-240V
Input Current	5A (RMS) @ 115V AC, 2.5A (RMS) @ 230V AC
Frequency	50-60Hz
Consumables	SEEKER Cartridges and LSD Reagent Kits
Software	Spot Logic*
Communication	USB 2.0
Data Output	CSV, XML, HTML

\*Assumes 2 SEEKER workstations (8 SEEKER instruments), 3 runs per day  
 (40 samples x 4 instruments per workstation x 2 workstations x 3 runs = 960 total samples)

CONTACT US

+1 (919) 891-0432 | [www.baebies.com](http://www.baebies.com)

© 2017 Baebies, Inc.